

## INTRODUCTIONS

### 01 01 TARGETED BODY SCULPTING

Lutronic introduces the AccuSculpt, the first 1444nm wavelength device for laser-assisted lipolysis. The system has received FDA regulatory clearance and utilizes a pulsed Nd:YAG laser to deliver a 1444nm wavelength to accurately and effectively remove fat. It is safe for all areas of the body and “targets fat cells more precisely producing less scattering and less collateral heat to surrounding tissue,” says Daniel Man, MD, a board-certified plastic surgeon in Boca Raton, Florida. “Fat is removed from fatty compartments without hurting other important structures like nerves and vessels. It significantly reduces post-operative discomfort...patients often return to work and other activities in a few days.” Contact: 888.588.7644, [www.lutronic.com](http://www.lutronic.com).

### 02 02 POSTPROCEDURE CARE

Intense Recovery Treatment from REVALÉSKIN offers 1.5% antioxidant-rich CoffeeBerry extract plus green tea, grape seed extract and skin conditioners to soothe and hydrate dry, damaged, postprocedure skin. Safe for all skin types, the nonirritating formulation works to plump, restore and protect skin following noninvasive procedures and is also effective for everyday rejuvenation. Contact: 866.738.2539, [www.revalaskin.com](http://www.revalaskin.com).

### 03 03 SMALL FOOTPRINT, BROAD INDICATIONS

The JOULE Aesthetic Laser Platform from Sciton offers multiple laser technologies through three delivery modes—arm, fiber and broadband light (BBL). The upgradeable, modular system can be used to provide laser-assisted lipolysis, deep and fractional skin resurfacing, hair removal and vascular treatments. Current modules include the Contour TRL 2940nm laser for skin peeling; Profractional and Profractional-XC 2940nm laser for skin resurfacing; the 1064nm ClearScan for hair removal; 1319nm ThermoScan for fine lines, wrinkles, acne and acne scarring; SkinTyte infrared for tissue firming; BBL broadband light for vascular lesions, hair removal and acne; ProLipo Plus with 1064nm and 1319nm wavelengths for laser-assisted lipolysis; and 1319nm Pro-V for endovenous ablation. Contact: 888.646.6999, [www.joulelaser.com](http://www.joulelaser.com).

### 04 04 STABLE VITAMIN C

Biopelle introduces Ascoderm, a new L-ascorbic acid-based line of topical treatments developed to reduce the signs of aging. The line includes three contouring serums, a cleanser, two moisturizers and an eye cream. The products—which can be mixed-and-matched with a patient's daily skincare regimen—offer a stable vitamin C formulation that penetrates the stratum corneum and creates a persistent reservoir of L-ascorbic acid in the skin. In addition to vitamin C, the line includes green tea extract, vitis vinifera extract and ginkgo biloba. Contact: 866.424.6735, [www.biopelle.com](http://www.biopelle.com).

### 05 05 SMOOTH OUT CELLULITE

The new UniForm handpiece from Alma Lasers has received FDA clearance for temporary reduction in the appearance of cellulite. It combines unipolar radiofrequency to deliver therapeutic heat into the dermal and subdermal tissue, plus mechanical massage to increase circulation in the subcutaneous tissue. The UniForm utilizes Alma's In-Motion technology to gradually heat tissue reducing pain and the risk of injury while contracting both collagen and subcutaneous tissue for a firm, smooth appearance. It is the fifth module available on the Accent XL. Contact: 866.414.2562, [www.almalasers.com](http://www.almalasers.com).



## INTRODUCTIONS



### 06 NEW INJECTABLE OPTIONS

If your injectable patients are trying to squeeze by with less than optimal volume based on multiple syringe cost concerns, two new Radiesse syringe volumes from BioForm Medical can help you address their budget constraints without sacrificing outcomes. The Volume Advantage 1.5cc syringe allows you to treat patients who require more volume with one single syringe. For patients with moderate wrinkles and folds, the Moderate Fill syringe holds just 0.8cc for less waste. Contact: 866.862.1211, [www.bioform.com](http://www.bioform.com).

### 07 GUIDED LIPOLYSIS

Osyris Medical introduces the LipoControl for faster, safer laser-assisted lipolysis procedures. LipoControl features an onscreen electronic mapping system that allows you to see the distribution of heat in real time. The 890nm diode laser system automatically adjusts to compensate for the movement of the cannula—the laser power increases as the speed of the cannula increases and decreases as the cannula slows—to prevent burning. The system requires only tumescent anesthesia and offers continuous or pulsed—from 0.1 to 25 seconds—emission of energy through its 600micron fiber. Contact: 972.720.0425, [www.osyrismedicalusa.com](http://www.osyrismedicalusa.com).

### 08 INTENSE HYDRATION

SkinMedica Hydrating Complex offers both high and low molecular weight hyaluronic acid. The higher weight stays above the skin to provide instant hydration; the lower weight penetrates into dry skin to rebuild its moisture reservoir. The formulation includes three natural sugar compounds to stimulate the skin's natural production of HA and a hyaluronidase inhibitor to prevent the breakdown of HA and elastin for prolonged hydration and firmer skin. Contact: 866.867.0110, [www.skinmedica.com](http://www.skinmedica.com).

### 09 FASTER TREATMENT TIMES

Cover larger hair removal areas in less time with the LightSheer HS handpiece from Lumenis. With a spot size of 22mm x 35mm, the LightSheer HS treats larger areas—including the legs and back—in only 15 minutes while achieving greater depth of penetration. In addition, the handpiece incorporates vacuum-assisted technology to reduce pain and improve efficacy. The vacuum draws the skin into the handpiece, removing competing chromophores and pulling hair follicles closer to the energy source. Like the original LightSheer ET handpiece, the HS offers pulse widths up to 400ms. Contact 877.586.3647, [www.aesthetic.lumenis.com](http://www.aesthetic.lumenis.com).

### 10 JUST FOR MEN

Mesoestetic Men offers four products developed specifically for male dermatologic concerns including irritation after shaving, excess sebum production, dryness and rough skin texture. The lightweight cosmeceutical line features nonirritating formulations with a high concentration of actives. Products include Eye Contour High Repair to eliminate eye bags, dark circles and wrinkles; Hydra Revitalizer Power moisturizer; Aftershave Inhibitor Effect to soothe skin after shaving and delay hair regrowth; and Abdo Attack, a micro-impulse patch that attacks abdominal fat. Contact: 888.498.3842, [www.mesoesteticusa.com](http://www.mesoesteticusa.com).

